



Criteria	Subject Device	Comparable Device	Comparable Device	Conclusion
Device Name	HCT Tamper Evident Syringe Cap	IMI Tamper Evident Syringe Cap	BBraun Tamper Evident Syringe Cap	Same
Packaging Configuration	10 Caps per Tray	10 Caps per Tray	10 Caps per Tray	Same
Labeling	Instructions on how to open syringe cap trays and attach to syringes.	Instructions on how to open syringe cap trays and attach to syringes,	Instructions on how to open syringe cap trays and attach to syringes,	Same
Sterility	Sterile EO (10^-6)	Sterile EO (10^-6)	Sterile EO (10^-6)	Same
Number of Uses	Single Use, Rx only	Single Use, Rx only	Single Use, Rx only	Same
Cap Material	Polypropylene Inner Cap	Polypropylene Inner Cap	Polypropylene Inner Cap	Same
Tamper Evident Design Feature	Cap Separation Upon Opening	Cap Separation Upon Opening	Cap Separation Upon Opening	Same



Product Data Sheet: HCT Tamper Evident Syringe Cap

NOTICE: The information in this document is the property of Yangzhou WEIDELI Trade Co., Ltd.. Unauthorized use is prohibited.

Manufacturer:	Yangzhou WEIDELI Trade Co., Ltd			
Description:	Tamper Evident Cap (with female, luer lock fitting)			
HCT Part Number(s):	HCT10R (Red) HCT10W (White) HCT10B (Blue) HCT10K (Black)	HCT100R-5B (Red) HCT100W-5B (White) HCT100B-5B (Blue) HCT100K-5B (Black)	HCT500R-5B (Red) HCT500W-5B (White) HCT500B-5B (Blue) HCT500K-5B (Black)	
Unit of Measure:	500 Units per Case	1000 Units per Case	1000 Units per Case	
Packaging Configuration:	Sterile Primary Packaging: 10 Units per Tray	Sterile Primary Packaging: 10 Trays per Bag	Sterile Primary Packaging: 500 Units per Bag	
	Secondary Primary Packaging: 50 Trays per Box	Secondary Primary Packaging: 10 Bags per Box	Secondary Primary Packaging: 2 Bags per Box	
Sterilization:	Ethylene Oxide (EO)			
Biocompatibility:	Biocompatibility testing has been conducted on materials used in the manufacturing of this device.			
Endotoxin:	Manufacturing lots of the device are routinely monitored for Bacterial Endotoxins.			
DEHP:	This product is made from materials that do not contain Di(2-ethylhexyl)phthalate (DEHP)			
Latex Statement:	This product does not contain natural rubber latex.			
Physical and Chemical Properties:	The physical and chemical properties of the plastic used manufacturing of this product have been tested.			
BSE/TSE Statement:				

BSE/ISE Statement:


Materials:	Item(s) Inside Luer Lock Cap Clear Outer Shell Cover Tray	Material(s) Used Polypropylene (PP) Polylac MABS Tyvek Polyethylene Terephthalate Glycol (PETG)
510(k) Number:	K121803	

European Conformity:

Additional Notes:



Specifications of the HCT Tamper Evident Syringe Cap



Specifications of the HCT Tamper Evident Syringe Cap

