


**Product Specification Sheet  
for  
Secure® EVA Legless Bags  
with  
Female Screw Connector**

**Product Features:**

- This product is not made with natural rubber latex or DEHP
- Ports and port tubes contain TOTM plasticized PVC
- EVA container
- Non-DEHP bag, ports and tubing
- Latex free
- Fill leg with female screw connector and male screw threads
- 250 ml to 4000 ml
- For use with filling devices requiring bags to have a female screw connector



Part #	Size (ml)	Description	Qty/Case	Case Price
66205	250	Secure® EVA Legless Bag – Female Screw Connector	50	
66210	500	Secure® EVA Legless Bag – Female Screw Connector	50	
66215	1000	Secure® EVA Legless Bag – Female Screw Connector	50	
66225	2000	Secure® EVA Legless Bag – Female Screw Connector	50	
66230	3000	Secure® EVA Legless Bag – Female Screw Connector	50	
66235	4000	Secure® EVA Legless Bag – Female Screw Connector	50	

Optional Accessories

Part #	Description	Qty/Case	Case Price
58997	Non-Reopening Clamp	50/bag	

(All Secure® EVA Legless Bags have an attached non-reopening clamp)

Manufactured by:

**THE METRIX COMPANY**  
4400 Chavenelle Road  
Dubuque, IA 52002  
563-556-8800

Distributed by: