


PMIFR100



Disposable, Single Use Device. Priming: 28.2ml (Approx), Length 100" (254cm) Universal Spike, upper pinch clamp, dial type flow controller, male spin luer lock.

PART#	PART DESCRIPTION	MATERIAL	LENGTH
1	SPIKE & DRIP CHAMBER	ABS+PVC	
2	PINCH CLAMP	ABS	
3	FLOW REGULATOR	ABS	
4	ROTATING MALE LUER LOCK	ABS	
5	NON DEHP PVC TUBING	PVC(NON DEHP)	92±2CM
6	NON DEHP PVC TUBING	PVC(NON DEHP)	143±2CM

