


## HCT21

Fluid Transfer Peristaltic Tubing Set  
(Compatible with Baxa/Baxter Products)



No.	Description	Material	Dimension	ID	OD	Wall	Hardness
1	Male Luer & Cap						
2	Non-DEHP PVC Tubing	PVC	140cm	4.7mm	6.4mm	0.85mm	80A
3	Connector (Big Hole)	ABS					
4	Silicone Tubing	Silicone	22.5cm	3mm	7mm	2mm	59A
5	Connector (Small Hole)	ABS					
6	NON-DEHP PVC Tubing	PVC	160cm	2.6mm	3.6mm	0.5mm	80A
7	Pinch Clamp	ABS					
8	Luer Holder	ABS					
9	Male Luer & Cap	SMMA & PP					
10	Silicone Pump Insert	ABS + Silicone	25.4cm				