Peristaltic Tubing Products from HCT

Proudly Providing Pharmacy Solutions for Over 30 years

(HCT)

(HCT)

.....

Peristaltic Pump Sets & Single-Use Systems (SUS)

- HCT products are compatible with top brands as listed below
- Talk to your sales rep about product programs and availability
- In-stock products will ship the same day if order is placed before 4 PM EST



	Order No.	Description	Qty
	HCT11	Fluid Transfer Peristaltic Tubing Set with Vented Spike / Luer Lock (Compatible with Baxa/Baxter Products)	10/CS
	НСТ21	Fluid Transfer Peristaltic Tubing Set Luer Lock / Luer Lock (Compatible with Baxa/Baxter Products)	10/CS
	НСТЗІА	Fluid Transfer Peristaltic Tubing Extension Set with ABS Luer Lock to ABS Luer Lock Low Priming Volume (Compatible with Baxa/Baxter Products)	10/CS
Q.	НСТ87А	Fluid Transfer Peristaltic Tubing Repeater Extension Set with ABS Luer Lock to ABS Female Luer Lock (Compatible with Baxa/Baxter Products)	50/CS
a series	НСТ29	Fluid Transfer Peristaltic Tubing Set Add-On with Vented Spike / Female Luer Lock (Compatible with Baxa/Baxter Products)	50/CS
	НСТ550	Fluid Transfer Peristaltic Tubing Set Add-On with Male Luer Lock to Male Luer Lock - 60" (Compatible with Baxa/Baxter Products)	50/CS
	НСТРА01	PharmAssist Fluid Transfer Peristaltic Tubing Set - Single with Spike / Luer Lock Low Priming Volume (Compatible with ICU PharmAssist)	10/CS
	НСТРА02	PharmAssist Fluid Transfer Peristaltic Tubing Set - Double with Spike / Luer Lock Low Priming Volume (Compatible with ICU PharmAssist)	10/CS
	НСТРАОЗ	PharmAssist Fluid Transfer Peristaltic Tubing Set - Triple with Spike / Luer Lock Low Priming Volume (Compatible with ICU PharmAssist)	10/CS




HCT Custom Peristaltic Pump Sets & Single-Use Systems (SUS)

HCT Can Meet Your Complex Configuration Needs

HCT is available for consultation to meet your design requirements.

Let Us Show You Our Customization Capabilities, Contact Your Local Sales Rep

Below is a custom peristaltic pump set & single-use system (SUS) that HCT created for a client.



No.	Park Description	Material	Remarks
1	3L EVA BAG	SANDED	30um
2	COEXTRUSION TUBE+MALE LUER	PVC(NON DEHP)&EVA+ABS	ID:6.5MM / OD:8.5MM / LENGTH:3CM LUER:5.0MM
З	FEMALE LUER (X2)	ABS	6.4MM
4	SILICONE TUBE (X2)	SILICONE	ID:6.4MM OD:9.5MM /LENGTH:6.4CM
5	CONNECTOR	PP	7.9MM x 6.4MM
6	SILICONE TUBE	SILICONE	ID:4.8MM OD:8.0MM /LENGTH:110CM
7	MALE LUER	PC	6.4MM
8	H CONNECTOR	ABS	BLUE
9	CONNECTOR	PP	7.9MM x 14.2MM
10	SILICONE TUBE	SILICONE	ID:12.7MM OD:19.05MM /LENGTH:10.16CM
11	CAPSULE FILTER1	707 (PP housing PES 0.22 Barbed hos	Filtration area : 570cm2 e ends)
12	SILICONE TUBE (X2)	SILICONE	ID:6.4MM OD:11.2MM LENGTH:7CM
13	CAPSULE FILTER2	L707 (PP housing PES 0.22 Barbed hos	Filtration area: 2041cm2 e ends)
14	FLANGE	PP	
15	SILICONE GASKET	SILICONE	
16	CLAMP	ABS	
17	END CAP	PP	
18	Metal Clamp (2X)	304	(0010.5-505R) 105
19	Metal Clamp (2X)	304	(0012.3-706R) 123
20-23	Metal Clamp (1X)	304	
24	Metal Clamp (3X)	304	(0015.7-706R) 157
25	Metal Clamp (1X)	304	(0014.2-706R) 142
26	PE BAG WITH TIE (2X)	PE+NYLONE	

Compatible with Brands Like:

- Bausch
- Cole-Parmer
- Flexicon
- Fisherbrand
- & More