






## QI Medical EnviroTest™

- Ideal for use in critical surface and gloved fingertip sampling to meet USP <797> requirements. Unique pedestal cap and hinged paddle
- Stores at room temperature with long shelf life without refrigeration
- · Log sheet & gummed labels included
- Talk to your sales rep about product programs and availability
- Orders will ship the same day if the order is placed before 4 PM EST





EnviroTest™ Growth Media Paddles for Surface, Air, or Glove Fingertip Sampling					
Order No.	Description	Shelf Life/Months	Qty		
ET1000	Double sided TSA (trypticase soy agar) growth med paddles formulated to inactivate many antimicrobic compounds		10 Testers/CS		
ET3000	Double sided MEA (malt extract agar) formulated to inactivate many antimicrobial compounds with added antibiotics to prevent bacterial overgrowth		10 Testers/CS		

GroMed™ Media Contact Plates/Settling Plates - Formulated with Tryptic Soy Agar (TSA) Media				
Order No.	Description	Shelf Life/Months	Qty	
GM0601	60mm contact plates for microbiological sampling a surface and fingertips. TSA and supplemented with Polysorbate 80 and Lecithin neutralizers		10/CS	
GM0901	90mm settling plates for microbiological sampling o surface and fingertips. TSA and supplemented with Polysorbate 80 and Lecithin neutralizers		10/CS	

GroMed™ Incubators - Lab Cell Culture Incubator for Pharmacy or Lab Use				
Order No.	Description	Oty		
QL10-140	Incubator 13"Wx11"Dx15"H, 3 shelves, dial thermometer	1/CS		
QL12-140	Large incubator 19"Wx13"Dx21"H, 3 shelves, dial thermometer	1/CS		



## GroMed™ Growth Media, Test Kits, Vials, & More

- Gummed labels included. Log sheets included for most products
- Talk to your sales rep about product programs and availability
- Orders will ship the same day if the order is placed before 4 PM EST



	Order No.	Description	Shelf Life/Months	Oty
	GM0030	TSB Growth media, 3mL of standard strength in glass ampules	24	20 Ampules/CS
	GM0100	TSB Growth media, 10mL tube, screw cap and needle access port	24	20 Tubes/CS
	GM0101	FTM Growth media, 10mL tube, screw cap and needle access port	12	20 Tubes/CS
	GM0200	TSB Growth media, 20mL vial, needle/spike septum	24	20 Vials/CS
CONT.	GM0201	FTM Growth media, 20mL vial, needle/spike septum	12	10 Vials/CS
DO MOD	GM0203	Fluid D, 20mL vial, needle/spike septum	12	20 Vials/CS
	GM0252	TSB Growth media, 25mL vial, 2X, needle septum	24	10 Vials/CS
	GM1400	TSB Growth media, 100mL bag with needle access port and spike port	24	10 Bags/CS
	GM1500	TSB Growth media, 100mL bag with Male Luer Lock port and spike port	24	10 Bags/CS





	Order No.	Description	Shelf Life/Months	Oty
	GM2000	TSB Growth media, 100mL vial with septum for needle or dispensing pin access	24	10 Vials/CS
	GM2001	FTM Growth media, 100mL vial with septum for needle or dispensing pin access	12	10 Vials/CS
Service Control of the Control of th	GM3000	Soybean-Casein Digest Media, approximately 3 grams of nonsterile powder in foil pouch, makes 100mL of TSB media	48	20 Pouches/CS
	GM4000	PDA Potato dextrose agar, 10mL tubes, test for fungal contamination	12	20 Tubes/CS
	GM5000	TSB Growth media, 500mL bag, 1X strength, breakaway covered needle access port and spike port	24	10 Bags/CS
	GM6005	5cc TSB Growth media in 12cc luer, lock syringe	24	10 Syringes/CS
	GM6105	5cc FTM (Fluid Thioglycollate) Growth Media in 12cc luer lock syringe	12	10 Syringes/CS
	GM7020	PATT™, Personal Aseptic Technique test kit, 20mL vials, 100mL partially filled minibags, test procedure and results log	24	5/Vials - 5 Bags
	GM7030	PATT2™, Personal Aseptic Technique test kit, 3mL ampules, 20mL vials, 100mL partially filled minibags,test procedure and results log	24	5 Ampules 5 Vials - 5 Bags
Garre Mec	EV0200	20mL sterile empty vial for use in aseptic technique tests. Not for patient use	36	20 Vials/CS
	EV2000	100mL sterile empty vial for use in aseptic technique tests. Not for patient use	36	10 Vials/CS



## **QI Medical Contamination Testers**

- Gummed labels included. Log sheets included for most products
- Talk to your sales rep about product programs and availability
- Orders will ship the same day if the order is placed before 4 PM EST





TuffTest2™ System - Aliquot Test for Emulsions and Suspensions				
Order No.	Description	Shelf Life/Months	Qty	
GM8100	TuffTest2™ contamination tester, 100mL TSB with TSA in sterile plastic bottle, needle access port	12	4 Testers/CS	
GM8101	TuffTest2™ contamination tester, 100ml FTM with Clostridial agar in sterile plastic bottle, needle access port	12	4 Testers/CS	

	QuickTest™ System - Full Filtration for All LVPs, Including Emulsions				
	Order No.	Description	Shelf Life/Months	Oty	
	QT1000	Contamination tester for use with:  • ADComp™ Compounder  • Automix Compounder (Clintec)  • Autocomp 6™ Compounder  • Exacta-Mix™ Compounder (Baxa)  • Nutrimix™ Compounder (Abbott)	24	10 Testers 10 100mL Bags of media	
	QT1100	Contamination tester with universal spike and 18" extension set	24	10 Testers 10 100mL Bags of media 10 16 ga. Needles	
E E	QT2000	Contamination tester for use with:  • Any pump set with male luer connector  • Repeater™ Pump (Baxa)  • Hyperformer™ Compounder (McGaw)  • PharmaAssist™ Pump (Excelsior)	24	10 Testers 10 100mL Bags of media	
	QT3000	Contamination tester for use with: Pinnacle compounder (BBraun)	24	10 Testers 10 100mL Bags of media	



## **QI Medical Testing Systems & More**

QT Micro™ Systems - Full Filtration for Small Volume Sterile Solutions				
Order No.	Description	Shelf Life/Months	Qty	
 TM6000	0.22µ Contamination tester, female luer lock inlet, female luer lock outlet, TSB Media	24	10 Testers 10 5cc Syringes of media	
TM6001	0.22µ Contamination tester, female luer lock inlet, female luer lock outlet, FTM Media	12	10 Testers 10 5cc Syringes of media	
TM6100	0.22µ Contamination tester, female luer lock inlet, male luer slip outlet, TSB Media	24	10 Testers 10 5cc Syringes of media	
TM6101	0.22µ Contamination tester, female luer lock inlet, male luer slip outlet, FTM Media	12	10 Testers 10 5cc Syringes of media	

ChemoTest™ - Hazardous Drug Handling Test Kits				
	Order No.	Description	Shelf Life/Months	Qty
Co Maded	СТ9530	Validate the technique of health care workers who manipulate hazardous drugs. Uses fluorescein dye and a UV light to detect breaks in safe technique and effectiveness of hazmat clean up procedures. 3mL and 50mL vials	24	5 Tests

PALL Syringe Filters					
Order No.	Description	Qty			
HP1002	0.22µ filter, 25mm, sterile, female luer lock inlet, male luer lock outlet, hydrophilic	50 Filters/CS			
PG005	Bubble Point Tester for use with syringe filters, 0-100psi range	1 Tester/CS			