

pfmmedical

*prime*MIDLINE[™]

Power injectable, low-profile Midline catheters for reliable peripheral access or selected intravenous therapies, blood sampling and power injection of contrast media.

Talk to your sales rep about product consultation, product programs and availability.



Features & Benefits

Low-Profile Luer Hubs:

- · Lays flat for patient comfort under occlusive dressings
- Allows for easier grasping and twisting when changing connections
- Resists cracking and deformation

3F High Performance Offering:

- Accommodates smaller vessels
- Achieve high flow performance at 5ml/sec
- Minimize thrombotic risks
- True to Size Catheter Length with No Taper at the Distal End

Soft Polyurethane Material:

- Provides ease of insertion
- Body softens to minimize risk of damage to the inner vessel wall

High Injection Rates Provide Optimal Imaging for Contrast Studies:

- 5 ml/sec for 3F Single Lumen Catheter
- 7 ml/sec for 4F and 5F Single and 5F Dual Lumen catheters

Order No.	Lumens	French Size	Catheter Length	Lumen Gauge	Max CT Rate (ml/sec)	Max Pressure (psi)	Qty
PFM3SML20	Single	ЗF	20cm	19G	5.0	300	5/CS
PFM4SML20	Single	4F	20cm	18G	7.0	300	5/CS
PFM5SML20	Single	5F	20cm	17G	7.0	300	5/CS
PFM5DML20	Dual	5F	20cm	18G/18G	7.0	300	5/CS

Kit Components

- Midline Catheter
- Tuohy Borst Side Port with Pre-Loaded Stylet
- Echogenic Introducer Needle 21G x 7cm
- Guidewire 45cm

- Tearaway Introducer
- Syringe 10cc
- Safety Scalpel
- Tape Measure
- Neutral Needleless Valve
- Statlock Plus Securement Device *Statlock is a registered trademark of C.R. Bard, Inc or an affiliate

<mark>pfm</mark>medical



Achieving the performance of power injection for short- and long-term peripheral access therapies up to 5 mL/sec. The primePICC catheters combine power injection and standard PICC therapies into one proficient catheter.

Talk to your sales rep about product consultation, product programs and availability.



Features & Benefits

- primePICC catheter with Tuohu Borst Adapter and Stylet (Stylet is not included in Over-the-Wire Kits)
- 21G Safety Introducer Needle with Echo Tip (21G Introducer Needle with Echo Tip in Over-the-Wire Kits)
- Neutral Displacement Needle-Free Connector (1 per Lumen)
- Statlock[®] Securement Device
- Tape Measure
- Premium Nitinol Guidewire with Floppy Tip
- 10 cc Syringe
- Tearaway Introducer

• Safety Scalpel

Available in 3 Configurations							
Order No.	Lumens	French Size	Lumen Gauge	Max Flow Rate (ml/sec)	Qty		
Basic PICC Kits							
PFM2CT3S	Single	ЗF	21G	5.0	5/CS		
PFM2CT4S	Single	4F	19G	5.0	5/CS		
PFM2CT5S	Single	5F	17G	5.0	5/CS		
PFM2CT4D	Dual	4F	20G/20G	4.0	5/CS		
PFM2CT5D	Dual	5F	19G/19G	5.0	5/CS		
PFM2CT5T	Triple	5F	22G/19G/22G	5.0	5/CS		

HCT)

Nursing & Over-the-Wire Configurations							
Order No.	Lumens	French Size	Lumen Gauge	Max Flow Rate (ml/sec)	Qty		
Nursing PICC Kits							
PFM2CT3SN	Single	ЗF	21G	5.0	5/CS		
PFM2CT4SN	Single	4F	19G	5.0	5/CS		
PFM2CT5SN	Single	5F	17G	5.0	5/CS		
PFM2CT4DN	Dual	4F	20G/20G	4.0	5/CS		
PFM2CT5DN	Dual	5F	19G/19G	5.0	5/CS		
PFM2CT5TN	Triple	5F	22G/19G/22G	5.0	5/CS		

Over-the-Wire PICC Kits

PFM2CT3SL	Single	ЗF	21G	5.0	5/CS
PFM2CT4SL	Single	4F	19G	5.0	5/CS
PFM2CT5SL	Single	5F	17G	5.0	5/CS
PFM2CT5DL	Dual	5F	19G/19G	5.0	5/CS
PFM2CT5TL	Triple	5F	22G/19G/22G	5.0	5/CS



