



AirLife Capnography
Cannulas



Proudly Providing [Airway](#)
Solutions for Over 30 years

- Reflective connector compatible with leading capnography monitors
- Reflective ETCO2 connector & thread grip O2 connector
- 40% Less dead space in facepiece
- Hydrophobic 0.22 micron filter



Item numbers with an “-L-” feature nafion tubing which is designed for moisture management; ideal for long term patient use.

An example of nafion tubing shown below.



EtCO2/O2 Filtered Capnography Cannula

Size: Adult
Interface: Oral/Nasal

Order No.	ETCO2 Length	O2 Length	Qty
4MSF1-7-13-25	13'	7'	25/CS
4MSF1-7-6-25	6'	7'	25/CS
4MSF1-L-7-13-25	13'	7'	25/CS

Features: Pigtail

4MSF1-0-13-25	13'	2"	25/CS
4MSF1-0-6-25	6'	2"	25/CS

Features: Nafion Tubing

4MSF1-L-7-6-25	6'	7'	25/CS
----------------	----	----	-------



EtCO2/O2 Filtered Divided Cannula

Size: Pediatric
Interface: Oral/Nasal

Order No.	ETCO2 Length	O2 Length	Qty
4MSF2-PED-7-6-25	6'	7'	25/CS

Features: Nafion Tubing

4MSF2-L-PED-7-13-25	13'	7'	25/CS
4MSF2-L-PED-7-6-25	6'	7'	25/CS



EtCO2 Filtered Cannula

Size: Adult
Interface: Nasal


Order No.	ETCO2 Length	O2 Length	Qty
4MSF3-6-25	6'	-	25/CS
Features: Nafion Tubing			
4MSF3-L-6-25	6'	-	25/CS



EtCO2/O2 Filtered Capnography Cannula

Size: Adult
Interface: Nasal

Order No.	ETCO2 Length	O2 Length	Qty
4MSF3-7-13-25	13'	7'	25/CS
4MSF3-7-6-25	6'	7'	25/CS
Features: Nafion Tubing			
4MSF3-L-7-6-25	6'	7'	25/CS



EtCO₂ Filtered Cannula

Size: Pediatric
Interface: Nasal

Order No.	ETCO ₂ Length	O ₂ Length	Qty
-----------	--------------------------	-----------------------	-----

Features: Nafion Tubing

4MSF3-L-PED-6-25	6'	-	25/CS
------------------	----	---	-------



EtCO₂/O₂ Filtered Capnography Cannula, Nafion Tubing

Size: Infant
Interface: Nasal

Order No.	ETCO ₂ Length	O ₂ Length	Qty
-----------	--------------------------	-----------------------	-----

4MSF3-L-INF-7-6-25	6'	7'	25/CS
--------------------	----	----	-------



EtCO₂/O₂ Filtered Capnography Cannula

Size: Pediatric
Interface: Nasal

Order No.	ETCO ₂ Length	O ₂ Length	Qty
-----------	--------------------------	-----------------------	-----

4MSF3-PED-7-6-25	6'	7'	25/CS
------------------	----	----	-------

Features: Nafion Tubing

4MSF3-L-PED-7-6-25	6'	7'	25/CS
--------------------	----	----	-------



EtCO₂ Filtered Capnography Sampling Line

Order No.	ETCO ₂ Length	O ₂ Length	Qty
-----------	--------------------------	-----------------------	-----

4MSF8-13-25	13'	-	25/CS
4MSF8-6-25	6'	-	25/CS

Features: Nafion Tubing

4MSF8-L-6-25	6'	-	25/CS
--------------	----	---	-------



EtCO₂ Filtered Sampling Line with Adapter

Order No.	ETCO ₂ Length	O ₂ Length	Qty
-----------	--------------------------	-----------------------	-----

4MSF5-13-25	13'	-	25/CS
4MSF5-6-25	6'	-	25/CS

Features: Nafion Tubing

4MSF5-L-13-25	13'	-	25/CS
4MSF5-L-6-25	6'	-	25/CS

Talk to your sales rep about product programs and availability.

