







Deliver medication quickly and painlessly with the DART. Atomized medication delivery results in therapeutic onset delivery similar to intravenous delivery. There's no special training required for a system wide adoption. The delivery method is painless for the patient, helping them remain calm in a high stress situation.

- Measure and deliver with syringe accuracy of (dart100): ± 5%
- Approximately 50 Microns
- NO IV line needed
- Available for all ages
- ALS and BLS approved







	DART <sup>™</sup> Intranasal Mucosal Atomization Device	
Order No.	Description	Qty
DART100	DART Nasal Device with 3mL Syringe	10/CS
DART110	DART Nasal Device with 1mL SyringePediatric	10/CS
DART140	DART Nasal Device with 3mL SyringePediatric, vial adapter	10/CS
DART300	DART Nasal Device without Syringe Neonatal	10/CS






Deliver medication quickly and effectively for difficult airways with the Reach. With the flexible stylet, adjust the tube for a customized approach.

## Ideal for:

- Airway Management; Difficult Airway, Awake Intubations
- Endoscopy
- Trans Esophageal Echocardiograms
- Bronchoscopy







DART-Reach Laryngo-Tracheal Mucosal Atomization Device		
Order No.	Description	Oty
REACH600 REACH700 REACH720	DART-Reach Device with 5mL Syringe DART-Reach-(218mm) Device without 5mL Syringe DART-Reach-(117mm) Device without 5mL Syringe	25/CS 25/CS 25/CS

- Delivers a smooth consistent atomized mist of medication
- · Flexible stylet gives you control over your procedure
- ISO Standard luer lock
- · No prepackaged medication to suit with your protocols