


Pulmodyne® **EZ-Spray™**

The EZ-Spray was designed for rapid anesthetization of patients. It is a single patient use, closed system, preventing the possibility of cross-contamination. Adjust the oxygen source to operate at 7-10 LPM according to desired flow rate. The 180° directional tip provides the ability to effortlessly direct and deliver medication to the upper pharynx. It is a quick and thorough anesthetization process.

- Disposable & single patient use
- Cross-contamination prevention
- Features a directional tip



EZ-Spray 313-100



EZ-Spray™ Powered Atomization Device with Directional Tip

Order No.	Description	Qty
313-100	EZ-Spray (60cc) with 7' Oxygen Tubing (shown above)	10/CS
313-101	EZ-Spray (15cc) with 7' Oxygen Tubing	10/CS
313-102	EZ-Spray SMR (15cc)	10/CS