


Latex Free, for Single Patient Use, & Disposable



· Allows clinician to deliver consistent PIP and PEEP pressures, which ensures greater safety when ventilating babies' fragile lungs, particularly in an emergency situation

• Integrated pressure manometer provides clear visualization of ventilatory pressure ranges

- Relief valve helps prevent over-pressurization of lungs
- · Able to provide CPAP
- Includes a cushioned face mask for proper sealing during mask ventilation
- · Easy to read controller and manometer makes set up simple and quick
- Universal red Fits-All O2 connector allows immediate visualization of proper oxygen connection



SafeT w/Pop-Off Valve		
Order No.	Size	Oty
ST710	X-Small Neonate	10/CS
ST711	Small Neonate	10/CS
ST712	Small Infant	10/CS
ST713	Infant	10/CS
ST714	Neonate Round	10/CS
ST715	Infant Round	10/CS
ST717	Without Mask	10/CS

